

Brushless Electric Wheelchair

This wheelchair is with a **smart designed frame**, the **Net weight is 23kg only**. It has brushless motor and controller, this not only minimize the weight of the chair, but also increase the driving distance.



YATTLI INDUSTRY CO., LTD.

Tel: +86 755 817 071 99

Web: www.yattll.com



Email: info@yattll.com

Factory: No.2 Lingchang Road, PingShan Village, TangXia Town, DongGuang City, GuangDong, China

Office: RM.1403, BLK.B of Duocai Office Complex, GuanlanHigh-tech Zone, Longhua DIST., Shenzhen, China

Parameter

Front wheel: 8" wheel

Rear Wheel: 12.5" wheel

Net Weight: 23 kg

Brake: Electromagnetic Braking system

Batteries: Lithium Battery 24V 11Ah

Motors (pair): Brushless Motor 180W

Charger: 24V 2A/3A/5A

Controller: Reliable Brushless Controller

Maximum Range: 15-18 km

Suspension System: Yes

Ground Clearance: 8 cm

Max Speed: 6 km/h

Charging Hour: 5-6 Hours

Max Slope: 12 degrees


Weight Capacity: 135 KG



Features

New Brushless type light weight power wheelchair; Weight only 23kg without battery and supports passengers weighing up to 135kg; New designed wheel with absorbing ability; Suspension system on back of frame offers super comfort; Ergonomics shaped seat, removable for easy cleaning.

Optional accessories



Folding step



Footnotes

- 1) Range and speed vary with userweightterrain type, battery charge battery condition and tyre pressure.
- 2) Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
- 3) Battery weight may vary based on manufacturer.

Color

Seat back pad color



Frame color

