

KIDS YE200C

— Electric Wheelchair —



PRODUCT PARAMETER

Model No.	YE200C
Seat width	39 cm
Seat depth	38.5 cm
Front wheel	8"PU solid wheel
Rear wheel	12.5"PU solid wheel
Battery	6Ah/9Ah *2pcs
Motor	Brushless 180W
Charger	24V 2A/3A/5A
Controller	Brushless controller
Brakes	Electromagnetic braking
Drive distance	15-18 km
Recline	Yes
Weight capacity	135 kg
Max speed	8 km/h
Net weight	32.8kg
Expanded Size	105*54*110 cm



YATTL INDUSTRY CO.,LTD



Tel: +86 755 81707199

Email: info@yattll.com

Website: www.yattll.com

Factory: YALU TECHNOLOGY PARK, NO. 88 Pingshan Industrial Rd, Tangxia, Dongguan City, Guangdong, China.

Office: Rm 907, Yichuang International Centre, No. 8 Gaoerfu Blvd, Longhua Dist, Shenzhen, Guangdong, China

YE200C Kids



Remote controlled reclining



Dual battery design



Safety protection baffle



Suspension system



Dual 180W brushless motors

Comfortable Reclining

Quick Released Battery



Product Details



8"PU solid wheel



12.5"PU solid wheel



Intelligent remote controller



Safety protection baffle

Features

- With multi-angle reclining backrest.
- Suspension system on the back frame.
- Waist-guard cushion offers better care to kids.
- Unique battery designed easily removable.
- Remote control for easier operating.