

## PRODUCT PARAMETER

Model No.	YE200R
Seat width	45 / 52 / 56 cm
Seat depth	43 cm
Front wheel	8" solid wheel
Rear wheel	12.5" Solid, Aluminum rim
Battery	6Ah/9Ah *2pcs
Motor	Brushless motors 250W*2
Charger	24V 2A/3A/5A
Controller	Brushless controller
Brakes	Electromagnetic braking
Drive distance	15-18 km
Recline	Yes
Weight capacity	150 kg
Max speed	6 km/h
Net weight	31.6kg
Expanded Size	110*60*98 cm

# TO BE A YATLL AGENT

## 2024

### Requirements

1. At least 1 lawful and registered local shop or sale place managed by yourself.
2. Legally registered company.
3. Have some knowledge and strong interest in electric wheelchairs

### Our Support

1. Agency area protection.
2. The best quality and competitive price.
3. Prior rights to the new products with new opportunity.
4. No minimum order quantity.
5. Soonest response on any of your questions.
6. Free parts supply policy.

JOIN US



YATLL INDUSTRY CO.,LTD



**Tel:** +86 755 81707199

**Email:** info@yattll.com

**Website:** www.yattll.com

**Factory:** YALU TECHNOLOGY PARK, NO. 88 Pingshan Industrial Rd, Tangxia, Dongguan City, Guangdong, China.

**Office:** Rm 907, Yichuang International Centre, No. 8 Gaoerfu blvd, Longhua Dist, Shenzhen, Guangdong, China

MODEL:YE200R

# Super Comforts Electric Wheelchair



Reclining



Remote control



Dual battery design

# Electric Recline

## Product Details



Recline by remote control



Angle-adjustable headrest



Extended footpad



Puncture proof wheel

Brushless controller  
with recline function



Extended footpad



105°

150°



12 Ah / 18Ah optional



Unique Dual-battery Design



Electric/Manual Mode switch

## Features

1. Electric recline & Wider seat & 150kg
2. Quick and easy removable battery, Lithium battery 6Ah X 2.
3. Brushless motor 24V 250W each.
4. Ergonomic-designed seat cushion with reflective stripes
5. Angle-adjusted headrest.
6. Extendable and comfortable footrest.
7. With Medical motor and Cylinder to support smooth Reclining of backrest.